Vision
■ ESTHETIC 2

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Introduction

Who is standing still is not getting better and will not move forward.
In order to counteract this, we the Wohlwend AG have moved and want to introduce The new Vision-Esthetic 2.
Retaining the advantages of Vision Esthetic, to facilitate processing was our goal. Reduce to the maximum.
This in this case is our motto.

What advantages does a “novelty” of a proven system for you?

- The processing (handling) is better
- A smaller range
- The alignment of all our Presssystemes in a color scheme

The lively color reproduction and natural fluorescence of Vision Esthetic 2 correspond to the natural teeth. In both ways, a simple, as well in an extended layering, with a wide range of Incisals and Transpa masses. There is a great flexibility in the full ceramic range. There are ingots with different opacities available which can be used for staining and layering techniques.

Compatibility with alloys in a CTE range 15,9- 16,7 (25-500°C)
Antagonists friendly Vickers hardness 470
Highest color accuracy at V shades A1-D4
Easy handling
High stability
Excellent aesthetics

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Material Properties

We confirm that all porcelain powders stated on this list have been tested and conform to EN ISO 9693:2000 and 6872:1998.





Keramik Ceramic Powder	Löslichkeit Solubility max. 100 µg/cm ² (+5)	Biegefestigkeit Flexural Strength min. 50 MPa	CTE (25/500°C) [ppm/K]		TG [°C] ± 10	Korngröße D90/10 [µm] Particle Size Distribution	Crystallinity detected by XRD
			2 Brand Firing	4 Brand Firing			
Opaker Opaker modifier	35	> 100	13,0	13,0	560	30	amorphous Glass/ Leucite
Dentin Opaque Dentin	16	> 75	14,2	14,4	480	60	amorphous Glass/ Leucite
Schneide Enamel	16	> 75	14,2	14,4	480	60	amorphous Glass/ Leucite
Schulter Shoulder	25	> 75	13,9	14,1	500	60	amorphous Glass/ Leucite
Clear, Transpa	16	> 75	14,2	14,4	480	60	amorphous Glass/ Leucite
Malfarben Stain	-	-	14,2	14,4	480	60	amorphous Glass/ Leucite
Glasur Glaze	-	-	14,2	14,4	480	60	amorphous Glass/ Leucite

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Preparation guide

- 1. Crowns:**
Reduce the incisal edge by 2 mm. Labial/Buccal areas 1–1,5 mm. Distinct chamfer or 90° shoulder. Palatinal minimum space of 1 mm. No sharp edges.
- 2. Inlays/Onlays:**
Inner surfaces should not be parallel, but open with an angle of 15° to the occlusal and approx. 2 mm deep. No feather edges nor opposing bite at tooth/inlay junction.
- 3. Veneers/Laminates:**
All preparation margins as chamfers of 0.5 mm. Minimal thickness 0.5 mm to the labial anterior and up to 1 mm to the incisal edge.

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

Note the alloy parameters

- Compatible with alloys having a CTE of 15,9-16,7(25–500°C).
- Ideally at 16,2(25–500°C).
- If the CTE is higher 16,4 (25–500°C):
Use slow cooling cycle!
- If below 16,1 (25–500°C):
Remove restoration from lift as soon as furnace has fully opened.
- These parameters may vary between different alloys.
- If you have second thoughts we recommend to use Vision-Bonding.

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Metal processing prior to porcelain application

1. Process alloys according to manufacturer's instructions.
 Avoid sharp edges. Unless specified otherwise, the alloy should be sandblasted with min. 110µm AL2O2 at 2 Bar.
 Make a oxid firing according to manufacturers instruction.
 A uniform oxide layer is desired.
2. In most cases, especially on NP, the oxides are to be removed using min. 110µm aluminous oxide.
 (Refer to manufacturers manual!)
 A clean and grease-free surface must be present.
3. Apply Vision-bonding according to the manual and fire it.
4. The result after the firing of the bonder.

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First step: 1st Opaque

1. Apply opaque with 70% semi-covering layer. It's the technician's own device to use a brush or a glass spatula.

After first firing. Shiny appearance after firing.

Attention: Refer to firing chart.
Changed firing cycle for NP alloys.

Second step: 2nd Opaque

2. Apply opaque with covering layer.

Egg shell appearance after firing.

Options

To achieve internal effects, characterizations can be made by using opaque modifiers, IO".

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**Third step:
 Shoulder**

The shoulder powders A, B, C, D and Neutral are giving you the chance to create a shoulder for every shade situation.

Here is the combination chart:

Shouldercombinationchart

Can change by using an individual colour scheme !!!

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Erster Dentinbrand.

Fourth step: Build up

1. In case of less space apply a thin layer of opaque-dentin (occlusal, cervical, etc.).
2. Build up with dentin (shade that is given) in full contour.
3. Make a small cut back from incisal and approximal. Apply the borders and mamelons with a thin layer of transpa neutral.
4. Complete the anatomic form of the tooth with enamel according to the shade. A light oversize of the build up, brings a perfect result, due to the small shrinkage during the firing cycle. Use a mixture of chroma-dentine and dentin (50/50) to fill up interdental space.

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Fifth step: Correction

1. If needed correct the form with Transpa Clear.
2. Proceed with second Dentine firing.

Sixth step: Finishing

When the contact points and the occlusion are worked out, you have got two possibilities to glaze. Glaze firing without glaze powder and subsequent a mechanical polish (special rubber wheels, brush, Dia glaze and so on). This is only recommended if the situation is well known.

Glaze firing with glaze powder. Mix it with glaze liquid and apply a thin layer. For individual characterization use Vision-Esthetic 2 stains.

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Vision-Esthetic 2 full-ceramic

Crowns:

- Reduce the incisal edge by 2 mm.
- Labial/Buccal areas 1–1.5 mm.
- Distinct chamfer or 90° chamfer.
- Palatinal minimum space of 1 mm.
- No sharp edges.

Inlays/Onlays:

Inner surfaces should not be parallel, but open with an angle of 15° to the occlusal and approx. 2 mm deep. No feather edges nor opposing bite at tooth / inlay junction.

Veneers/Laminates:


All preparation margins as chamfers of 0.5 mm. Minimal thickness 0.5 mm to the labial anterior and up to 1 mm to the incisal edge.

Staining technique:

For inlays, onlays, veneers and single crowns. Restorations that have been waxed to full contour and pressed using transparent ingots are completed by using the appropriate stain.

Layering technique:

The layering technique for maximized esthetic values. Pressed copings are characterized using the appropriate dentine and incisal overlay powders including transparents and opal may also be used to obtain an even higher level of esthetics.

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Wax up

- Use only for full ceramic recommended 100% ash free waxes and sprues. Do not use debubbilizer!
- Apply spacer in a thin layer 1 mm above the prep margin on crowns, veneers, inlays and onlays.
- Staining technique: Full contour wax up. No loss of dimensions during pressing and firing processes.
- Build up technique: Either full wax up, reduced wax up or waxed as a coping. Minimum thickness 0.8 mm or 2/3 press material with 1/3 overlay ceramic.
- Spruing: Use 6 or 8 gauge ash free sprues with a length of 4–7 mm. Coping are sprued with an outward taper at both sprue base and coping. Sprues are connected to the coping in a straight line so as not to restrict the flow of the molten ingot. Sprues should be angled slightly outward with a 3 mm distance between copings.

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Investing

- Wax weight up to 0.6 g = 1 ingot
- Wax weight 0.6–1.3 g = 2 ingots

Optional:

Paper ring former or silicon ring former. Remove strip from paper, form a ring to the appropriate mark. Slip onto base and adapt stabilizing ring. (These steps are not necessary when using silicon ring). Pour special investment to just below the stabilizing ring. (Slightly over the mark in the silicon ring). Carefully remove the stabilizing ring and fit the ring former, twisting it carefully. Investment must leak from the opening.

After bench set, remove paper or silicon ring, base and ring former, then clean the bottom surface with a plaster knife. This must be absolutely flat, so that the ring stands perfectly upright in the pressing furnace.

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x¹ Layering technique, only in case of suitable preparation and sufficient space conditions!

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



Devesting

- Mark the length using spare plunger. Cut along the line, using seperatingdisk.
- Attention, use mask and dust suction.
- Carefully break segments apart.
- Use glass beads only. Start with 60 psi pressure.
- As soon as objects appear turn pressure down to 35 psi proceed blasting carefully.
- Do not use aluminium oxide.
- Carefully separate the sprue from the objects with a sharp diamond disk.

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Finishing

Carefully separate the sprue from the objects with a sharp diamond disk. Check the fitting and finish with fine diamonds burs and silicon rubber wheels. Do not use carbide burs or high speed hand pieces without water cooling.

Staining technique

Prepare objects to desired contour (texture, contacts, etc.). Sandblast with 50 μ aluminium oxide with low pressure and clean.

Additions can be made, using overlay material. Customize the restoration according to situation, using stains. This can be repeated several times.

After desired effects are achieved, glaze once or twice as required. Firing is carried out on Vision firing paste.

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Build up technique

Thoroughly clean the reduced, possibly individualized copings. Complete the reduced or missing portions, using the equivalent overlay materials. A correction cycle can be carried out if necessary. Prepare objects to the final contour (texture, contacts, etc.). Now glaze the restoration and use stains to add small nuances if desired. If you don't want to work with glaze powder, the material is also easy to polish mechanically.

Respect minimum thickness of 0.8 mm!

Fire on firing pillow or corresponding pins, molars are only to be fired on firing pillow.

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Firing chart PFM

Recommended firing parameters.
 Temperatures can vary due to different furnace performances.

* Final temperature for NP: 1562°F

If alloy CTE is higher than 16,4 (25–500°C) use a slow cooling cycle of 3 minutes.

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Firing chart full-ceramic

Recommended firing parameters.

Temperatures can vary due to different furnace performances.

Staining / Internal staining technique:

Schichttechnik:

Schichttechnik Layering technique	Starttemperatur Start temperature	Trockenzeit Dry time	Temperaturanstieg Heat rate	Endtemperatur End temperature	Haltezeit Hold	Vakuum Vacuum	Erscheinungsbild Appearance
Dentin 1 Dentine 1	450°C 842°F	6 min.	45°C/min. 113°F/min.	770°C 1418°F	1 min.	ja yes	leicht glänzend slightly shiny
Dentin 2 Dentine 2	450°C 842°F	6 min.	45°C/min. 113°F/min.	760°C 1400°F	1 min.	ja yes	leicht glänzend slightly shiny
Glanz ohne Glasurmasse Glaze without glaze powder	450°C 842°F	4 min.	45°C/min. 113°F/min.	770°C 1418°F	1 min.	–	glänzend shiny
Glanz mit Glasurmasse Glaze with glaze powder	450°C 842°F	4 min.	45°C/min. 113°F/min.	750°C 1382°F	1 min.	–	glänzend shiny

Layering technique:

Sortimente

TESTKIT NR.1 (klein)

O-A3 20g; D-A3 20g; OD-Honey 20g; I-59 20g; TR-Clear 20g;
Opaquer Liquid 25ml; Modelling Liquid 25ml;

TESTKIT NR.2 (groß)

5 Presspellets CT3; 5 Presspellets A3;
D-A3 20g; I-59 20g; TR-Clear 20g; Glaze 20g;
Stain : A 3g; White 3g; Redbrown 3g
1 Modelling Liquid 25ml; 1 Opaquer Liquid 25ml;
1 Glaze&Stain Liquid 25ml;
5X100gr Investment, 50ml Liquid

SORTIMENT S1

10 Opaquer Pulver 20g:
A1;A2;A3;A3.5; B1;B2;B3; C1;C2; D2
1 Opaquer Liquid 25ml

SORTIMENT S2

10 Dentin 20g : A1;A2;A3;A3.5; B1;B2;B3; C1;C2; D2
1 Modelling Liquid 25ml

Assortments

TESTKIT NR.1 (small)

O-A3 20g; D-A3 20g; OD-Honey 20g; I-59 20g; TR-Clear 20g;
Opaquer Liquid 25ml; Modelling Liquid 25ml;

TESTKIT NR.2 (large)

5 Pressables CT3; 5 Pressables A3;
D-A3 20g; I-59 20g; TR-Clear 20g; Glaze 20g;
Stain : A 3g; White 3g; Redbrown 3g
1 Modelling Liquid 25ml; 1 Opaquer Liquid 25ml;
1 Glaze&Stain Liquid 25ml;
5X100gr Investment, 50ml Liquid

ASSORTMENT S1

10 Opaquer powder 20g:
A1;A2;A3;A3.5; B1;B2;B3; C1;C2; D2
1 Opaquer Liquid 25ml

ASSORTMENT S2

10 Dentin 20g : A1;A2;A3;A3.5; B1;B2;B3; C1;C2; D2
1 Modelling Liquid 25ml

SORTIMENT S3

4 Incisal 20g : 57;58;59;60
7 Transpa 20g : Clear;Neutral;Opal;White;Red;Yellow;Grey
1 Glaze 20g
1 Modelling Liquid 25ml; 1 Glaze Liquid 25ml

SORTIMENT S4

4 Int.Opaquer 20g : White; Orange; Yellow; Brown;
2 Opaque-Dentin 20g : White;Honey
4 Chroma Dentin 20g : A;B;C;D
1 Gingiva Darkpink 20g
5 Shoulder powder 20g : A;B;C;D,Neutral
1 Opaquer Liquid 25ml; 1 Modelling Liquid 25ml
1 Shoulder Liquid 25ml

MALFARBEN-SORTIMENT

A 3g; B 3g; C 3g; D 3g;
White 3g; Honey 3g; Oliv 3g; Redbrown 3g; Chestnut 3g;
Grey 3g; Blue 3g; Glaze 3g,
1 Glaze&Stain Liquid, 1 Pinsel 00

ASSORTMENT S3

4 Incisal 20g : 57;58;59;60
7 Transpa 20g : Clear;Neutral;Opal;White;Red;Yellow;Grey
1 Glaze 20g
1 Modelling Liquid 25ml; 1 Glaze Liquid 25ml

ASSORTMENT S4

4 Int.Opaquer 20g : White; Orange; Yellow; Brown;
2 Opaque-Dentin 20g : White;Honey
4 Chroma Dentin 20g : A;B;C;D
1 Gingiva Darkpink 20g
5 Shoulder powder 20g : A;B;C;D,Neutral
1 Opaquer Liquid 25ml; 1 Modelling Liquid 25ml
1 Shoulder Liquid 25ml

STAIN-ASSORTMENT

A 3g; B 3g; C 3g; D 3g;
White 3g; Honey 3g; Oliv 3g; Redbrown 3g; Chestnut 3g;
Grey 3g; Blue 3g; Glaze 3g,
1 Glaze&Stain Liquid, 1 Brush 00

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Accessory

- Porcelain croco
- Brushes
- Mixing trays
- Ceramic isolation pen
- Firing paste
- Plunger investment
- Plunger Aluminium oxide
- Silicon ring former
- Investment
- Bonding

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WOHLWEND AG
Wohlwend AG Dental Manufaktur