Vision
■ LOW

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TECHNICAL MANUAL

WOHLWEND AG

Einführung

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Introduction

Competence means for us to take responsibility for reaching better solutions in a complete system. How do we get an upgrade in color scheme in its reproduction? Standard procedure to increase the efficiency. The requirement for the technician today is: better results with less complexity. So we want to hit this target with you.

Wohlwend AG, striving to meet the demands of the world market has expanded its ceramic line with Vision Low. Built on the advantages of Vision Esthetic, Vision Low offers greater flexibility as it is compatible with a wide range of alloys in a standard CTE range. Dual Low fusing porcelain for pressable ceramic and PFM restorations.

- Compatible with PFM alloys CTE range 13.8–14.9 (25–500°C)
- Press to metal CTE range 13.8–14.2 (25–500°C)
- Kinder to natural dentition Vickers hardness 470
- Accurate shade reproduction with V shades A1–D4
- Easy handling
- No slumping during build up
- Outstanding esthetics

The vital shade reproduction and natural fluorescence of Vision Low correspond to natural teeth. This replication of nature can be obtained by either using a basic two layer build-up or by incorporating additional incisal and transpas into the restoration. And nowadays a great flexibility in full ceramics. Ingots with different opacity for a usage in Stains-Layering and Press over Metal technique.

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Material Properties

We confirm that all porcelain powders stated on this list have been tested and conform to EN ISO 9693:2000 and 6872:1998.





Keramik Ceramic Powder	Löslichkeit Solubility max. 100 µg/cm ² (± 5)	Biegefestigkeit Flexural Strength min. 50 MPa	WAK CTE (25/500°C)		TG ± 10	Korngröße D90/10 Particle Size Distribution	Crystallinity de- tected by XRD
			2 Brand Firing	4 Brand Firing			
Opaquer, Opaquer Modifier	35 µg/cm ²	> 100 MPa	13.0 ppm/K	13,0 ppm/K	560°C	30 µm	amorphous glass / leucite
Dentine, Opaque Dentine	16 µg/cm ²	> 75 MPa	11.4 ppm/K	11.8 ppm/K	510°C	60 µm	amorphous glass / leucite
Schneide Enamel	16 µg/cm ²	> 75 MPa	11.4 ppm/K	11.8 ppm/K	510°C	60 µm	amorphous glass / leucite
Schulter Shoulder	25 µg/cm ²	> 75 MPa	12.2 ppm/K	12.2 ppm/K	530°C	60 µm	amorphous glass / leucite
Clear, Transpa	16 µg/cm ²	> 75 MPa	11.4 ppm/K	11.8 ppm/K	510°C	60 µm	amorphous glass / leucite
Malfarben Stain	–	–	11.4 ppm/K	11.8 ppm/K	520°C	20 µm	amorphous glass / leucite
Glasuren Glaze	–	–	11.4 ppm/K	11.8 ppm/K	520°C	20 µm	amorphous glass / leucite

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Preparation guide

1. Crowns:
Reduce the incisal edge by 2 mm. Labial/Buccal areas 1–1,5 mm. Distinct chamfer or 90° shoulder. Palatinal minimum space of 1 mm. No sharp edges.
2. Inlays/Onlays:
Inner surfaces should not be parallel, but open with an angle of 15° to the occlusal and approx. 2 mm deep. No feather edges nor opposing bite at tooth/inlay junction.
3. Veneers/Laminates:
All preparation margins as chamfers of 0.5 mm. Minimal thickness 0.5 mm to the labial anterior and up to 1 mm to the incisal edge.

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
Alloy Parameters

- Compatible with alloys having a CTE of 13.8–14.9 (25–500°C).
- Ideally at 14.2 (25–500°C)
- If the CTE is higher 14.5 (25–500°C):
Use slow cooling cycle!
- If below 14.1 (25–500°C):
Remove restoration from lift as soon as furnace has fully opened.
- These parameters may vary between different alloys.
- If you have second thoughts we recommend to use Vision Bonding.

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Metal processing prior to porcelain application

1. Process alloys according to manufacturers instructions. Avoid sharp edges. Unless specified otherwise, the alloy should be sandblasted with min. 110 µm aluminum oxide at 29 psi.
2. A uniform oxide layer is desired.
3. In cases with NP alloys, the oxides are to be removed using min. 110 µm aluminum oxide. (Refer to manufacturers manual!). A clean and grease-free surface must be achieved. Now apply Vision Bonding according to the manual.

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First step: 1st Opaque

1. Apply opaque with 70% semi-covering layer. It's the technicians own device to use a brush or a glass spatula.
2. After first firing. Shiny appearance after firing.


Attention: Refer to firing chart.
Changed firing cycle for NP alloys.

Second step: 2nd Opaque


3. Apply opaque with covering layer.
4. Egg shell appearance after firing.

Options

To achieve internal effects, characterizations can be made by using opaque-modifiers „IO“.

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
Third step: Shoulder

1. Three different shoulder powders are giving you the chance to create a shoulder for every shade situation.
2. To lower the values mix in incisal bleach.
3. To lower the chromas mix in transpa clear.

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





Fourth step: Build up

1. In case of less space apply a thin layer of opacious dentin (occlusal, cervical, etc.).
2. Build up with dentin (shade that is given) in full contour.
3. Do a small cut back from incisal and approximal. Apply the borders and mamelons with a thin layer of transpa neutral.
4. Complete the anatomic form of the tooth with enamel according to the shade. A light oversize of the build up, brings a perfect result, due to the small shrinkage during the firing cycle.

Use some opaque dentine to fill up interdental space. Fire first Dentine firing.

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Fifth step: Correction

1. After first dentin firing.
2. If needed correct the form with Transpa Clear.
3. Proceed with second Dentine firing.

Sixth step: Finishing

4. When the contact points and the occlusion are worked out, you have got two possibilities to glaze.

Glaze firing without glaze powder and subsequent a mechanical polish (special rubber wheels, brush, Dia glaze and so on). This is only recommended if the situation is well known.

Glaze firing with glaze powder. Mix it with glaze liquide and apply a thin layer. For individual characterization use Vision Low stains.

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Vision Low Full-ceramic

Crowns:

- Reduce the incisal edge by 2 mm.
- Labial/Buccal areas 1–1.5 mm.
- Distinct chamfer or 90° chamfer.
- Palatinal minimum space of 1 mm.
- No sharp edges.

Inlays/Onlays:

Inner surfaces should not be parallel, but open with an angle of 15° to the occlusal and approx. 2 mm deep. No feather edges nor opposing bite at tooth / inlay junction.

Veneers/Laminates:

All preparation margins as chamfers of 0.5 mm. Minimal thickness 0.5 mm to the labial anterior and up to 1 mm to the incisal edge.

Staining technique:

For inlays, onlays, veneers and single crowns. Restorations that have been waxed to full contour and pressed using transparent ingots are completed by using the appropriate stain.

Layering technique:

The layering technique for maximized esthetic values. Pressed copings are characterized using the appropriate dentine and incisal overlay powders including transparents and opal may also be used to obtain an even higher level of esthetics.

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Wax up

- Use only for full ceramic recommended 100% ash free waxes and sprues. Do not use debubbilizer!
- Apply spacer in a thin layer 1 mm above the prep margin on crowns, veneers, inlays and onlays.
- Staining technique: Full contour wax up. No loss of dimensions during pressing and firing processes.
- Build up technique: Either full wax up, reduced wax up or waxed as a coping. Minimum thickness 0.8 mm or 2/3 press material with 1/3 overlay ceramic.
- Spruing: Use 6 or 8 gauge ash free sprues with a length of 4–7 mm. Coping are sprued with an outward taper at both sprue base and coping. Sprues are connected to the coping in a straight line so as not to restrict the flow of the molten ingot. Sprues should be angled slightly outward with a 3 mm distance between copings.

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Investing

- Wax weight up to 0.6 g = 1 ingot
- Wax weight 0.6–1.3 g = 2 ingots



Optional:

Paper ring former or silicon ring former.

Remove strip from paper, form a ring to the appropriate mark. Slip onto base and adapt stabilizing ring. (These steps are not necessary when using silicon ring). Pour special investment to just below the stabilizing ring. (Slightly over the mark in the silicon ring).

Carefully remove the stabilizing ring and fit the ring former, twisting it carefully. Investment must leak from the opening. After bench set, remove paper or silicon ring, base and ring former, then clean the bottom surface with a plaster knife. This must be absolutely flat, so that the ring stands perfectly upright in the pressing furnace.

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Layering technique

(Optional for staining technique)

x¹ Layering technique, only in case of suitable preparation and sufficient space conditions!

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
Approximate value for calibrated press furnaces

Start-temp.	Heizrate Heat rate	Endtemp. Final temp.	Haltezeit Hold time	Presszeit Press time	Vakuum Vacuum	Druck Pressure
°C °F	°C/min. °F/min.	°C °F	min.	min.	%	Bar psi
700 1292	60 140	950 1742	20	10	100	4,5-5 65-72

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



Devasting

- Mark the length using spare plunger. Cut along the line, using separating disk.
- Attention, use mask and dust suction.
- Carefully break segments apart.
- Use glass beads only.
Start with 60 psi pressure.
- As soon as objects appear turn pressure down to 35 psi proceed blasting carefully.
- Do not use aluminium oxide.
- Carefully separate the sprue from the objects with a sharp diamond disk.

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Finishing

Carefully separate the sprue from the objects with a sharp diamond disk. Check the fitting and finish with fine diamonds burs and silicon rubber wheels. Do not use carbide burs or high speed hand-pieces.

Staining technique

Prepare objects to desired contour (texture, contacts, etc.). Sandblast with 50 μ aluminium oxide with low pressure and clean. Additions can be made, using overlay material. Customize the restoration according to situation, using stains. This can be repeated several times. After desired effects are achieved, glaze once or twice as required. Firing is carried out on Vision firing paste.

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
Build up technique

Thoroughly clean the reduced, possibly individualized copings. Complete the reduced or missing portions, using the equivalent overlay materials. A correction cycle can be carried out if necessary. Prepare objects to the final contour (texture, contacts, etc.). Now glaze the restoration and use stains to add small nuances if desired. If you don't want to work with glaze powder, the material is also easy to polish mechanically.

Respect minimum thickness of 0.8 mm!
Fire on firing pillow or thin metal pins, molars are only to be fired on firing pillow.

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Firing chart PFM

Recommended firing parameters. Temperatures can vary due to different furnace performances.

* Final temperature for NP: 1562°F

If alloy CTE is higher than 14.5 (25–500°C) use a slow cooling cycle of 3 minutes.

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Firing chart full-ceramic





Recommended firing parameters.
Temperatures can vary due to different furnace performances.

Staining / Internal staining technique:

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Press over Metal

Follow the alloy parameters.

The CTE range of the alloys should be in-between 13,8–14,2 (25–500°C). Prepare the surface, refer to the alloy manufacturer's instruction.

Opaque

Use powder or paste opaque from the Vision Classic ceramics.

1st Opaque

Apply opaque with (70%) semi-covering layer. Bake with 970°C. Shiny appearance after firing.

2nd Opaque

Apply opaque with covering layer. Bake with 950°C. Egg shell appearance after firing.

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Wax up

Use 100% ash free waxes only like in staining technique.

Full contour wax up. No loss of dimension during pressing.

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Pressing

Use transparent dentin pellets for press over metal technique. Devesting and finishing as usual.

Staining technique

Individual nuances can now be added with Vision Low stains.

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16 Dentine powder A1–D4 20 g
2 Modellingliquid 25 ml

ASSORTMENT S4

Incisal: 57, 58, 59, 60, 20 g
Opalincisal: Neutral, Blue 20 g; Intensivincisal: White, Yellow 20 g
Transpa: Clear, Neutral, Red, Yellow, Grey, Glaze 20 g,
1 Modellingliquid 25 ml,
1 Glaze & Stainliquid 25 ml

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ASSORTMENT S6

Shouldermasses: Creme, Yellow, Redbrown 20 g, 1 Shoulderliquid 25 ml

ASSORTMENT S7

Pressable ingots A1-D4, D-BA00, D-BBO, D-BBOO, 5 of each

ASSORTMENT S8

Pressable ingots CL, NT, NW, T1, T2, T3, T4, CT1, CT2, CT3, CT4, OP1, OP2, 5 of each

SORTIMENT S9

Pressable ingots TA1, TA2, TA3, TB1, TB2, TB3, TC1, TC2, TD2, Bleach, 5 of each

SORTIMENT S10

Bleach Dentine: BL1, BL2, BL3, Incisal-Bleach, Modellingliquid 25 ml
 Bleach opaque, Opaqueliquid 25 ml

STAIN-ASSORTMENT

White, Honey, Olive, Redbrown, Chestnut, Grey, Blue, Glaze, A, B, C, D, 3 g
 1 Glaze & Stainliquid 25 ml

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Accessory

- Porcelain croco
- Brushes
- Mixing trays
- Ceramic isolation pen
- Firing paste
- Plunger investment
- Plunger Aluminium oxide
- Silicon ring former
- Investment
- Bonding

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