
Vision  
■ ZIRKON

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TECHNICAL MANUAL

WOHLWEND AG

## Einführung

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## Introduction

Innovations afford new materials! With its high strength and delicate frameworks, zirconium oxide is enjoying increased popularity. Wohlwend's Vision Zirkon provides a new standard for shade accuracy. Its exact CTE settings and excellent adhesion result in optimum performance on different zirconium oxide frameworks. This protects restorations from cracks, checks and fractures.

With our new press to zircon system we once more unite high aesthetics with constant quality and easily work engineering. The new development is compact and applicable for every dental technician. Valuable time is saved through the singular in fully contour principle and less needed for the rework. Hereby we achieve a higher efficiency with a high-quality result.

- Shade accuracy
- Superb natural esthetics
- Easy handling
- Great stacking / carving abilities
- Outstanding adhesion and excellent bonding strength
- Exact CTE settings CTE 9.5

Press over Zirkon

- pressable with or without ceramic shoulder
- natural color depth and translucency
- Stain color system in 4 V-Shades
- only one press liner
- Cut Back technique with Vision Zirkon possible

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## Material Properties

We confirm that all porcelain powders stated on this list have been tested and conform to EN ISO 9693:2000 and 6872:1998

Klasse Class	Löslichkeit Solubility max. 100 µg/cm <sup>2</sup> (± 5)	Biegefestigkeit Flexural Strength min. 50 MPa	WAK CTE (25/500°C)		TG ± 10	Korngröße D90/10 Particle Size Distribution	Crystallinity de- tected by XRD
			2 Brand   Firing	4 Brand   Firing			
Liner	16 µg/cm <sup>2</sup>	> 70 MPa	9,5 ppm/K	9,5 ppm/K	550 °C	60 µm	amorphous Glass
Dentin Opak-Dentin	16 µg/cm <sup>2</sup>	> 70 MPa	9,5 ppm/K	9,5 ppm/K	550 °C	60 µm	amorphous Glass
Incisal Int.Incisal	16 µg/cm <sup>2</sup>	> 70 MPa	9,5 ppm/K	9,5 ppm/K	550 °C	60 µm	amorphous Glass
Shoulder	16 µg/cm <sup>2</sup>	> 70 MPa	9,5 ppm/K	9,5 ppm/K	580 °C	45 µm	amorphous Glass
Transparent	16 µg/cm <sup>2</sup>	> 70 MPa	9,5 ppm/K	9,5 ppm/K	550 °C	60 µm	amorphous Glass
Stain	–	–	9,5 ppm/K	9,5 ppm/K	550 °C	20 µm	amorphous Glass
Glaze	–	–	9,4 ppm/K	9,4 ppm/K	550 °C	20 µm	amorphous Glass

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

To reach perfect results it is required to prepare the framework exactly and according to the manufacturers instructions.

### First step: Liner

A thin layer liner according to the shade will be put on the framework and fired, before starting second step.

### Second step: Shoulder

Three different shoulder powders are giving you the chance ...

- to create a shoulder for every shade situation.
- to lower the values mix in intensive Incisal white.
- to lower the chroma mix in transpa clear.

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

Erster Dentinbrand.

### Third step: Build up

1. In case of less space apply a thin layer of opacious dentin (occlusal, cervical, etc.).
2. Build up with dentin (shade that is given) in full contour.
3. Make a small cut back from incisal and approximal. Apply the borders and mamelons with a thin layer of transpa neutral.
4. Complete the anatomic form of the tooth with enamel according to the shade. A light oversize of the build up, brings a perfect result, because of the small shrinkage during the firing cycle.

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### Fourth step: Correction

1. After first dentin firing.
2. Work out the crown and correct the form with Transpa Clear.

To achieve opalescence in the incisal area, TR Opal can be mixed into all transpa materials.

Even on a simple build up, natural light transmission and reflection can be achieved.

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
## Fifth step: Finishing

1. When the contact points and the occlusion are worked out, you have got two possibilities to glaze.
2. Glaze firing without glaze powder and a subsequent mechanical polish (special rubber wheels, brush, Dia glaze and so on). This is only recommended if the situation is well known.


Glaze firing with glaze powder. Mix it with glaze liquid and apply a thin layer.

For individual characterization use Vision Zircon stains.

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## Press over Zircon

It is possible to press over all commercially Zirconia systems (green, sintered or hiped). The grinded or sintered Zirconia framework has to be prepared with special grinders or rubberwheels to manufacturer's specifications. Here, the manufacturer's minimum strengths of copings and thickness of connections have to be followed.

### First step: Liner

1. Apply the special universal liner in a thin and watery consistence to the framework and fire as shown in the firing chart. This serves to reduce the brightness value of the zircon framework material, and to get a better bonding.

**Attention:** Different firing temperature.

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
## Second step: Wax up

Use 100% ash free waxes only! Put the framework without pressure on the isolated model. It is possible to wax up in full contour, or to reach a higher esthetic work in a reduced way (cut back).

In such cases apply from our Vision Zircon incisals and transpa powders. There is no loss of dimension during pressing.

## Third step: Spruing

Use 6 or 8 gauge ash free sprues with a length of 4 to max. 7 mm. Sprues are connected to the coping in a straight line so as not to restrict the flow of the molten ingot. Sprues should be angled slightly outward with a 3 mm distance between copings.

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#### Wichtig:

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### Fourth step: Investing

Wax weight up to 0.6 g = 1 Ingot  
Wax weight 0.6–1.3 g = 2 Ingots

#### Optional:

Paper ring former or silicon ring former.

Remove strip from paper, form a ring to the appropriate mark. Slip onto base and adapt stabilizing ring. Pour special investment to just below the stabilizing ring.

Carefully remove the stabilizing ring and fit the ring former, twisting it carefully. Investment must leak from the opening. After bench set, remove paper / Silicon ring, base and ring former, then clean the bottom surface with a plaster knife. This must be absolutely flat, so that the ring stands perfectly upright in the pressing furnace.

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## Fifth step: Pressing


For full contour press technique use Transpa Dentin pellets. For cut back technique use Dentin pellets.




Look for pressparameters and combinationcharts.

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## Sixth step: Devasting

Mark the length using spare plunger. Cut along the line, using separating disk. Attention, use mask and dust suction. Break segments apart.

Carefully devaste the objects by sandblasting. Use glass beads only.

Start with 60 psi pressure, as soon as objects appear turn pressure down to 35 psi proceed blasting carefully. DO NOT USE ALUMINIUM OXIDE!

Carefully separate the sprue from the objects with a sharp diamond disk.

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## Neunter Schritt: Glaser

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## Seventh step: Layering (cut back)

After working out the press channels, layer up with Vision Zircon. Complete the incisal part and fire according to the chart. With dentin you can layer supplemented portions. If necessary make a correction firing. Before glazing you can customize using Vision Zircon stains. Fire on ceramic pins. It is important to ensure that the ceramic shoulder has no contact with the pin. It can be glazed or polished.

## Eighth step: Staining technique

Individual nuances can now be added with Vision Zircon Stains.

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## Zircon firing chart

Recommended firing parameters.

Temperatures can vary due to different furnace performances.

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16 Dentine powder A1–D4 20 g  
2 Modellingliquid 25 ml

### **ASSORTMENT S2**

Incisal: 57, 58, 59, 60, 20 g  
Opalincisal: Neutral, Blue 20 g; Intesiveincisal: White, Yellow 20 g  
Transpa Clear 20 g, Glaze 20 g, Liner: 1 to 3 and Pressliner 20 g,  
1 Modellingliquid 25 ml, 1 Glaze & Stainliquid 25 ml

### **ASSORTMENT S3**

Transpa: Neutral, Yellow, Blue, Grey 20 g  
Opaquedentine: White, Honey 20 g,  
Gingival: Pink and Softpink 20 g  
Shouldermasses: Creme, Yellow, Redbrown 20 g  
1 Modellingliquid 25ml, 1 Shoulderliquid 25 ml

### **ASSORTMENT S4**

Pressables each 5 pc: A1, A2, A3, B1, B2, B3, C1, C2, D2, Bleach,  
TA1, TA2, TA3, TB1, TB2, TB3, TC1, TC2, TD2

### **STAIN-ASSORTMENT**

White, Honey, Olive, Redbrown, Chestnut, Grey, Blue 3 g,  
Glaze, A, B, C, D, 3 g, Glaze & Stainliquid 25 ml

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## Accessory

- Porcelain croco
- Brushes
- Mixing trays
- Ceramic isolation pen
- Firing paste
- Plunger investment
- Plunger Aluminium oxide
- Silicon ring former
- Investment

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**Wohlwend AG